

102 SALINE IMPLANT SETUP		updated March 25, 2013
SEE CASE NOS.	107 Mastopexy with Saline Implants 108 Breast Augmentation with Saline Implants 109 BREAST IMPLANT REMOVAL and REPLACEMENT	
ROOM SETUP	Double Basin Set [1-2] liter bag Normal Saline, hang at top of bed sterile cysto tubing on field [2] blue one-way valve connectors [1] Twist-n-Ject stop cock [8-10] 60cc syringes [2] L irrigation saline on filed to fill sizers Kefzol irrigation, use blue bulb syringe to rinse sizers and implants	
IMPLANTS	Implant sizers - NOTE: as of 2013, sterile, non-reusable implant sizers will be provided by the OR. Implants - open only when requested, verify with Dr.C before opening. DO NOT REMOVE ANY STICKERS. KEEP ALL LABELS AND PACKING IN BOX. MARK BOX WITH <u>R</u> or <u>L</u> SIDE AND RECORD VOLUMES. INCLUDE <u>6</u> PT ID STICKERS FOR DR.C TO USE.	
SOLUTIONS IRRIGATION	! Normal saline irrigation [3] liters ! [2] Asepto syringes ! Kefzol irrigation for implants w/ blue bulb OPEN ONLY IF USING SALINE: ! Hang 2L sterile IV saline to fill implants ! Cysto tubing, single	
NEEDLES SYRINGES	[8-10] 60cc syringes for filling implants	
STERILE SUPPLIES SYRINGES and CONNECTORS	! [2] Asepto syringes ! [1] Blue bulb syringe for antibiotic irrigation OPEN ONLY IF USING SALINE: ! [8-10] 60cc syringes for filling sizers and implants ! [1] Twist-n-Ject stop cock ! [2] blue one-way valve connectors ! [1-2] Cysto tubing, single	
SETS and INSTRUMENTS TO OPEN	Double Metal Basin Set	
COMMENTS		
SEE ALSO, SAME AS	THIS IS A MASTER CARD FOR 107, 108, 113	
DATE MODIFIED	March 2006, March 2008, October 2010	