

[Back to Master List](#)

<h2 style="text-align: center;">Breast Reconstruction with Tissue Expander - First Stage</h2> <p>unilateral or bilateral, may include contralateral breast procedure such as breast reduction, augmentation, mastopexy, implant, mastectomy.</p>	
CATEGORY	Breast
CASE NO.	112 Breast Reconstruction with Tissue Expander - First Stage
OPERATIVE TIME	2-4 hours
USUAL LOCATION(S)	Pres Main
FREQUENCY	5/yr
ASSISTANT	Yes, Cheryl
ICD9 CODE:	174.1malignant Neoplasm Breast, 738.3 Acquired Deformity Chest
CPT CODE(S):	19357 INSERTION TISSUE EXPANDER
<p style="text-align: center;">ANESTHESIA NOTES</p> <p style="text-align: center;">PATIENT POSITION SETUP</p> <p style="text-align: center;">TABLE POSITION</p> <p style="text-align: center;">ROOM SETUP</p>	<p>GA, ETT or LMA, Kefzol 1gm IV</p> <p>supine, 6 chux, arms out on boards, covered with green towels and taped. Sequential teds, tattoo kit: 22g needle, kelly clamp, 3cc syringe 18g needle, 4x4, gloves</p> <p>pillow under knees, table in semi-recumbent position</p> <p>Pre-warm room and table</p>
<p style="text-align: center;">EQUIPMENT</p> <p style="text-align: center;">SPECIALS</p> <p style="text-align: center;">IMPLANTS</p> <p style="text-align: center;">PREP</p> <p style="text-align: center;">DRAPING</p> <p style="text-align: center;">CAUTERY</p> <p style="text-align: center;">SUCTION</p>	<p>Baer hugger, metric digital scale, small suction canisters for weighing tissue, bag decanter. Headlight and light source</p> <p>need fluid warmer for Gwen</p> <p>Implants - Call Jim Medina for implants. Open only when requested, verify with Dr. before opening. DO NOT REMOVE ANY STICKERS. KEEP ALL LABELS AND PACKING IN BOX. MARK BOX WITH SIDE AND RECORD VOLUME. INCLUDE 5 PT ID STICKERS FOR DR. TO USE.</p> <p>CHG prep set, down to blue chux</p> <p>Lap chole pack, mag pad, [3] quarter sheets, [2] mayo stand covers, [4-6] sticky paper drapes</p> <p>Bovie 50/50 pure with silver tip with Saf-T-Vac, Bovie extender</p> <p>[2] suction: regular with plastic Yankauer and Saf-T-Vac [2] suction tubing</p>
<p style="text-align: center;">SOLUTIONS IRRIGATION</p> <p style="text-align: center;">NEEDLES SYRINGES</p> <p style="text-align: center;">LOCAL ANESTHETIC</p> <p style="text-align: center;">MEDICATIONS</p>	<p>Normal saline [3] liters, [2] Asepto for normal irrigation Kefzol irrigation with blue bulb for implants</p> <p>[2] 20cc syringes, [2] long Quincke 22g spinal needles</p> <p>Solution "C" 200cc</p> <p><input type="checkbox"/> Bacitracin ointment 2 oz <input type="checkbox"/> methylene blue <input type="checkbox"/> Kefzol 1gm for irrigation</p>
<p style="text-align: center;">SETS TO OPEN</p> <p style="text-align: center;">INSTRUMENTS USED (ON MAYO STAND)</p> <p style="text-align: center;">■ FORCEPS</p> <p style="text-align: center;">■ SCISSORS</p>	<p>Plastic set, breast retractors. Lap set</p> <p>01 Adson, 03 Adson, 05 Adson, long Adson-Brown, Debakey</p> <p>Metzenbaum, Str Mayo, Curved Mayo</p>

	<p>■ NEEDLE HOLDERS</p> <p>■ RETRACTORS</p> <p>■ KNIFE BLADES</p> <p>■ OTHER</p>	<p>[2] Crile-Wood, [1] Webster, [1] Mayo-Hagar</p> <p>Facelift small and large, Army Navy, Richardson small and large, Deaver, breast retractors</p> <p>[2] No. 10, [2] No. 15</p> <p>Tonsil clamp, hemostat</p>				
	<p>STERILE SUPPLIES</p> <p>SYRINGES and CONNECTORS</p> <p>SPONGES</p> <p>NON-STERILE SUPPLIES</p>	<p>[2] marking pens with rulers</p> <p>[2] Asepto syringes, [2] non-disposable light handles</p> <p>[2] Asepto syringes, [1] Blue bulb syringe</p> <p>[1] Twist-n-Ject stop cock (leave all caps ON)</p> <p>[2] blue one way valve connectors</p> <p>[1] Cysto tubing</p> <p>[6] Lap sponges packs</p> <p>Hang 1L bag NS to fill expander</p>				
	SUTURES	●	No	Type	SN	status
		●	[2-4]	0-Vicryl CT-1		hold
		●	[2-4]	2-0 Vicryl CT-1	J945	open
		●	[2-4]	3-0 Vicryl SH	J416	open
		●	[2-4]	4-0 Vicryl PS-2	J496	open
		●	[2-4]	4-0 Monocryl PS-2	Y426	open
		●	[1-2]	3-0 Nylon PS-1	1663	open
		●	[]			
	DRAINS	[1-2] 15 Fr. Round Blake drains with <i>connectors and reservoir bulbs</i> - 2188.				
	DRESSINGS	<p>[2] tongue blades [2] ABD</p> <p>[2] large adaptic [2] 6" Ace</p> <p>[2] 4x4 [2] kerlix</p> <p>[2] cut drain sponges [2] safety pins</p>				
	POSTOPERATIVE	<p>Call for bed from 7D</p> <p>Lift with wrapping board, 3 people</p> <p>Wrapping board, large or small, baby powder</p> <p>Safety pin reservoirs to ace, mark 1 & 2 with Sharpie</p>				
	NOTES	May have specimen for pathology				
	<p>COMMENTS</p> <p>DIAGRAM</p> <p>DESCRIPTION</p> <p>SEE ALSO, SAME AS</p> <p>DATE MODIFIED</p>	<p>This operation is done after mastectomy for breast cancer. Tissue expanders are placed in the first stage and then the expansion is done as an outpatient. The second stage is done 3-4 months later in which the expander is removed and permanent saline implants are placed.</p> <p>April 30, 2006</p>				