

<b>204 REPAIR UNILATERAL or BILATERAL CLEFT LIP</b> <b>with or without REPAIR ANTERIOR HARD PALATE (AHP)</b> usually includes 220 PE TUBES BY DR. HORN, <u>always</u> includes 221 VIDEO SCOPE AND TOWER SETUP <span style="float: right;">updated June 5, 2016</span>	
ICD10 CODES: Q37.9 UNILATERAL CLEFT LIP & PALATE Q37.8 BILATERAL CLEFT LIP & PALATE	CPT CODES: 40700 REPAIR UNILATERAL CLEFT LIP 40701 REPAIR BILATERAL CLEFT LIP 42235 REPAIR ANTERIOR HARD PALATE
<b>ANESTHESIA NOTES</b>	<ul style="list-style-type: none"> <li>■ regular ET tube (no RAE)</li> <li>■ long circuit</li> <li>■ straight anesthesia connector</li> <li>■ Dr. C to start IV and sew in tube with 2-0 silk</li> <li>■ IV pump in room and charging</li> <li>■ temp probe</li> <li>■ Kefzol IV 25mg/kg; Decadron 0.5mg/kg</li> </ul>
<b>PATIENT POSITION   SETUP</b>	<ul style="list-style-type: none"> <li>■ Will turn bed counter clockwise after PE tubes</li> <li>■ Dr. C to suture tube with 2-0 silk SH,</li> <li>■ Rolled terry cloth towel under shoulders</li> <li>■ Large round gel doughnut under head (no cover)</li> <li>■ Pedi-wrap elbow restraints with D-ring (neonate size)</li> <li>■ Pedi Baer hugger #555</li> <li>■ blanket, green towels ■ 1" tape to secure tube</li> </ul>
<b>TABLE POSITION</b>	Lower foot, close head gap, DO NOT REMOVE HEAD PIECE Large round gel donut head rest. Place Baer after turning.
<b>ROOM SETUP</b>	Warm room 72°, warm bed and head rest w/ Baer Separate mayo stand for video equip
<b>EQUIPMENT and SPECIALS</b>	<ul style="list-style-type: none"> <li>● Dr. C light source w/headlight</li> <li>● Baer hugger</li> <li>● IV pump plugged in</li> <li>● Video tower and 2<sup>nd</sup> slave monitor</li> <li>● Omni Lapotrac scope holder and bed clamp</li> <li>● 3.5mm Endoscope, camera, light cord.</li> </ul>
<b>PREP</b>	PCMX with minimal water, cotton tip applicators to clean nose, wipe gently inside mouth, lightly over eyelids
<b>DRAPING and PACKS</b>	Cleft Palate Custom Pack, Towel Pack, White towels Head sheet (two green towels and a ¾ sheet), ¾ sheet, magnetic pad, rigid box for bovie, towel clips, ½" steri-strip for eyes, 1000 drape, 2" Kling for pharyngeal pack
<b>CAUTERY</b>	Bovie 8/8 pure, [1] Colorado tip, [1] 0118A Megadyne needle point tips, [2] holsters
<b>SUCTION x 2</b>	[1] 7fr Frazier tip <u>and</u> [1] Sontec peds fluted suction [1] 10fr NG suction catheter for end of case to suction stomach
<b>SOLUTIONS   IRRIGATION</b>	Kefzol irrigation 500mg/500cc NS, saline, water
<b>NEEDLES   SYRINGES</b>	[2] 3cc syringe, [1] 30cc syringe, [4] 30g needle 5/8", [1] 12g angiocath
<b>LOCAL ANESTHETIC</b>	<b>peds mixture:</b> 50/50 - Mix in 3cc syringe with 30g needle x 2 <ul style="list-style-type: none"> <li>• 1% lidocaine w/ epi 1:100,000 and</li> <li>• 0.5% Ropivacaine plain</li> </ul>
<b>MEDICATIONS</b>	<ul style="list-style-type: none"> <li>■ methylene blue in prep cup for tattooing lip</li> <li>■ Kefzol for irrigation</li> <li>■ Bacitracin ointment</li> </ul> <p style="text-align: center;">***for bleeding have available from anesthesia*** neosynephrine 0.2cc in plastic cup</p>

**204 REPAIR UNILATERAL or BILATERAL CLEFT LIP**  
*with or without REPAIR ANTERIOR HARD PALATE (AHP)*

page 2

**SETS TO OPEN** Dr. C Cleft Lip / Fine Instrument Set / Basin Set

**INSTRUMENTS USED (ON MAYO STAND) - 2 mayo stands**

■ <b>FORCEPS</b>	Packing forceps, [3] Cushings, .05 Adson, .5 TC gold CV fine, 0.5 silver Sontec CV fine, 0.5 titanium
■ <b>SCISSORS</b>	Gradle, fine, straight tenotomy, curved tenotomy, Iris, Stevens
■ <b>NEEDLE HOLDERS</b>	[4] Angled palate, [3] Castroviejos, smooth Webster
■ <b>RETRACTORS</b>	Millard-Dingman mouth gag, 2mm, 5mm, 10mm skin hooks, Ragnell
■ <b>KNIFE   BLADES</b>	[4] 30° super sharp blades #7531, [2] 15 BP scalpel blades
■ <b>ELEVATORS</b>	McCullough, Hockey stick, Tebbetts
■ <b>OTHER</b>	calipers, blue silicone cup, plastic med glasses

<b>STERILE SUPPLIES</b>	<ul style="list-style-type: none"> <li>■ [2] marking pens with ruler</li> <li>■ bacitracin ointment</li> <li>■ non-disposable light handles</li> <li>■ methylene blue</li> <li>■ med glasses</li> <li>■ defogger</li> <li>■ ½" steri-strips for eyes</li> <li>■ 10fr suction catheter</li> <li>■ [3] 3" cotton tip applic</li> <li>■ 7½ Protegrity latex gloves</li> </ul>
	HOLD ■ Gelfoam (surgifoam #1972) HOLD ■ Surgicel 2"x14" #1951
<b>SPONGES, etc.</b>	4x4 Raytec, 2" Kling for pharyngeal pack, [3] 3" cotton tip applicators, towel pack
<b>NON-STERILE SUPPLIES</b>	Pedi Baer hugger #555, Pedi-Wrap arm restraints neonatal

	No	Type	SN	status	
<b>SUTURES</b>	● [ 1 ]	2-0 Silk SH	K 833	open	
	● [ 8 ]	4-0 Chromic P3	1654	open 6	hold 2
	● [ 2 ]	5-0 Chromic P-3	687	open	
	● [ 2 ]	6-0 Chromic G-1single	796	open	
	● [ 2 ]	4-0 Monocryl RB1	Y214	hold 2	
	● [ 1 ]	5-0 Monocryl TF	Y433	open	
	● [ 1 ]	6-0 Monocryl P-1	Y489	open	
	● [ 1 ]	6-0 Prolene P-1	8697	hold	
	● [ 1 ]	7-0 PDS BV-1	Z 135	hold	
	● [ 2 ]	7-0 Prolene P-6	8648	open 1	hold 1

<b>DRESSINGS</b>	none - bacitracin to wound, remember to remove pharyngeal pack and cut suture, 10fr suction catheter to suction stomach.
------------------	--

<b>POSTOPERATIVE</b>	call for pediatric crib (6 <sup>th</sup> floor), 2 warm blankets in corral shape, connect IV infuser <u>before</u> moving baby, oxygen for transport
----------------------	--

<b>NOTES</b>	may not include anterior hard palate, may not have PE tubes, all cases will have video monitoring.
--------------	--

DATES MODIFIED	2006, 2007, 2008, 2009, Jan 2010, Aug 2012, Jan 2013, Sep 2015, Jun 2016
----------------	--